



Use Your Autosampler to Make Calibration Standards

Merlin K. L. Bicking, Ph. D.
ACCTA, Inc.


INTRODUCTION AND GENERAL SETUP

- We describe a feature of many autosamplers that produces dilutions and mixtures.
 - This approach is useful when:
 - High-value standard solutions must be prepared in a minimum volume with maximum accuracy.
 - Classic analytical dilutions involve large volumes of solvent and multiple volumetric flasks.
- Solutions Required
 - Stock Solution: Must be in compatible solvent
 - Receiving Vial(s):
 - Clean, empty, no cap
 - Vial inserts can be used for small volume dilutions
- HPLC System
 - Disconnect column
 - Use union or place column inlet tubing into waste container
 - Program pump to deliver desired dilution solution
 - Set flow to 0.5 - 1.0 mL/min
 - Any flow can be used, but the system pressure should be as low as possible.
 - Equilibrate system with flow through autosampler loop
 - "Mainpass" in Agilent systems

HOW DOES IT WORK?



EXAMPLE INJECTOR PROGRAM




OpenLab ChemStation C.01.09 - 1:100 Dilution

Needed to reset sampling syringe to home position.
This is the volume of stock solution.
Wash options will vary depending on your model.
This example uses a position relative to the sample. Depending on your model, you may have to specify a relative position for tray, row, and column.
Wash options will vary depending on your model.
Needed to complete the process.

- Other Method Settings
 - Pump
 - Flow: 1 mL/min
 - Stop Time: 0.1 min.
 - Injection
 - Settings for needle wash
 - Column
 - None needed
 - Detector
 - Default settings
- Injection/Sequence Settings
 - Set the sample as the location of the Stock solution
 - Or, you can set a fixed location in the program and the sample location is ignored.

CASE STUDY 1: CONSUMER PRODUCT INGREDIENT

- Method
 - Prepare 4 dilutions of a single stock solution in methanol for analysis of different products
 - 7%, 3%, 1%, and 5%
 - Pump Flow = 1 mL/min methanol
 - Locations
 - Stock Solution in P1-A1 (S)
 - 7% in P1-A3
 - 3% in P1-A4
 - 1% in P1-A5
 - 5% in P1-A6



Parameters for Case Study 1:

- 7%
- 3%
- 1%
- 5%


Results

Compare Replicate Preparations		Compare Preparation Levels		
	Avg Area	Level, %	Avg Area	RF
Prep 1	1574.65	7.00	1574.008	224.9
Prep 2	1578.40	3.00	670.017	223.3
Prep 3	1581.96	1.00	222.654	222.7
Avg	1578.3	5.00	1115.356	223.1
SD	3.65			223.0
RSD	0.23			0.345
			RSD	0.15

Compare with manually prepared standard

	Avg Area	Prep Time	Solvent Used
Manual Prep	1577.650	20 min	400 mL
Autosampler Prep	1578.337	12 min	15 mL
% Difference	0.04%		

CASE STUDY 2: DNPH DERIVATIVE



Autosampler (Program)		Manual (Serial)	
Lvl	Amt	Area	RF
1	0.1	17.231	172.31
2	0.2	34.431	172.16
3	0.5	85.098	170.20
4	1	171.45	171.45
5	2	343.7	171.85
6	5	859.47	171.89
7	10	1724.8	172.48
8	20	3449.2	172.46
9	50	8359.5	167.19
	Avg	171.3	
	SD	1.71	
	RSD	1.00	

Stock Vol./Total Vol. from Source

- 1000
- 50/1000 from 1000
- 20/1000 from 1000
- 10/1000 from 1000
- 5/1000 from 1000
- 1/1000 from 1000
- 50/500 from 10
- 0.5/500 from 5
- 0.2/500 from 2
- 0.1/500 from 10


Graphs showing Average Area vs Concentration for Autosampler and Manual methods.

Autosampler: $y = 173.28x + 85.173$, $R^2 = 0.9981$

Manual: $y = 167.57x + 18.285$, $R^2 = 0.9998$

CASE STUDY 3: MAKING MIXTURES

Prepare a mixture of three components, each diluted 1:25 from three individual stock solutions.



Sample location in Sequence Table is A1


Location is relative to sample, one column over (A2)

Location is relative to sample, one row down (B1)

Pump flow rate is 0.50 mL/min., so total volume is 500 µL

The Sample+ option allows you to specify a position relative to the Sample location. You specify the change for Tray, Plate, Row, and Column.

CASE STUDY 4: WIDE RANGE CURVE FOR GC



Level	Diluton	Conc	Avg Area	RF
1	10000	0.7992	368.600047	461.2
2	500	15.984	7275.203026	455.2
3	200	39.96	18492.79934	462.8
4	20	399.6	195377.5766	488.9
5	10	799.2	393765.4735	492.7
6	5	1598.4	803435.2665	502.6
7	2	3996	1993227.733	498.8
8	Stock	7992	4073564.138	509.7
			Avg	484.0
			SD	21.1
			RSD	4.4

Graph showing Average Area vs Concentration with regression equation: $y = 508.65x - 7977.8$, $R^2 = 0.9999$

SUMMARY

Limitations

- Small dilution factors (e.g., 1:2) are limited if large volumes are required. Vial inserts would be required e.g., 50 / 100
- UHPLC systems often have smaller injection limits. Requires use of inserts
- Some UHPLC pumps may not be as reliable at the low pressures observed with this method
- Mixing/agitation is required, so complete automation is not possible (yet). Future modifications may solve this problem.

Advantages

- Programmable autosamplers can be used to prepare mixtures and/or calibration standards with high accuracy
- Options will vary across models and software versions, but the general procedure is the same.
- Solvent use is significantly reduced
- Use of volumetric glassware is significantly reduced
- Only preparation of stock solution is required
- Smaller amounts of standards are required
- No need for storage of dilutions.
- Make what you need when you need it.
- Completion time is usually less than the corresponding manual method and does not require staff time during the program (in most cases).